


Product Datasheet

Rat MIP-1Alpha/CCL3 ELISA Kit (orb546944)

Description	Rat MIP-1Alpha/CCL3 ELISA Kit
Reactivity	Rat
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A24-A92
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate)
MW	10335 MW
Uniprot ID	P50229
Sensitivity	1 pg/ml
Expiration Date	Please enquire.



Rat CCL3 PicoKine ELISA Kit standard cur...