


## Product Datasheet

### Human CD93/C1qR ELISA Kit (orb546952)

<b>Description</b>	Human CD93/C1qR ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	78 pg/ml - 5,000 pg/ml
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Expression system for standard: NS0; Immunogen sequence: A24-K580
<b>Storage</b>	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
<b>Note</b>	For research use only
<b>Sample Types</b>	cell culture supernatants, serum, plasma ( heparin, EDTA) , saliva and urine.
<b>MW</b>	68560 MW
<b>Uniprot ID</b>	<b>Q9NPY3</b>
<b>Sensitivity</b>	10 pg/ml
<b>Expiration Date</b>	Please enquire.



Human CD93/C1qR PicoKine ELISA Kit stand...