


Product Datasheet

beta Tubulin Rabbit Monoclonal Antibody, HRP Conjugated (orb546964)

Description	beta Tubulin Rabbit Monoclonal Antibody, HRP Conjugated
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthesized peptide derived from human beta Tubulin
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	WB 1:5000-1:20000
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	50433 MW
Uniprot ID	Q13509
Sensitivity	10pg/ml
Expiration Date	12 months from date of receipt.



WB analysis of beta Tubulin expression i...