

Product Datasheet

GFP Rabbit Monoclonal Antibody, HRP Conjugated (orb546965)

Description	GFP Rabbit Monoclonal Antibody, HRP Conjugated
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthesized peptide
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	WB 1:3000-1:10000
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	25993 MW
Uniprot ID	Q96319
Expiration Date	12 months from date of receipt.

