


Product Datasheet


LDL Receptor Rabbit Monoclonal Antibody (orb547032)

Description	LDL Receptor Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthesized peptide derived from human LDL Receptor
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	WB 1:1000-1:2000IHC 1:100-1:500
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	95376 MW
Uniprot ID	P01130
Sensitivity	12pg/ml
Expiration Date	12 months from date of receipt.



WB analysis of LDLR expression in HepG2

...



Immunohistochemical analysis of paraffin...