



Product Datasheet

CD38 Rabbit Monoclonal Antibody (orb547137)

Description	CD38 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthesized peptide derived from human CD38
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	WB 1:500-1:2000IHC 1:100-1:500
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	34328 MW
Uniprot ID	P28907
Expiration Date	12 months from date of receipt.



WB analysis of CD38 expression in HeLa C...



Immunohistochemical analysis of paraffin...