



Product Datasheet

C Reactive Protein CRP Rabbit Monoclonal Antibody (orb547190)

Description	C Reactive Protein CRP Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	A synthesized peptide derived from human C Reactive Protein
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	WB 1:1000-1:2000IHC 1:500-1:1000IP 1:50
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	25039 MW
Uniprot ID	P02741
Expiration Date	12 months from date of receipt.



WB analysis of CRP expression in Etoposi...



Immunohistochemical analysis of paraffin...