


Product Datasheet

IP10 CXCL10 Rabbit Monoclonal Antibody (orb547207)

| | |
|----------------------------|---|
| Description | IP10 CXCL10 Rabbit Monoclonal Antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | A synthesized peptide derived from human IP10 |
| Form/Appearance | Liquid |
| Concentration | Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure. |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Application notes | WB 1:500-1:2000 |
| Isotype | Rabbit IgG |
| Clonality | Monoclonal |
| MW | 10881 MW |
| Uniprot ID | P02778 |
| Expiration Date | 12 months from date of receipt. |



WB analysis of IP10 expression in IP10 r...