



Product Datasheet

iNOS NOS2 Rabbit Monoclonal Antibody (orb547276)

Description	iNOS NOS2 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, IP, WB
Immunogen	A synthesized peptide derived from human iNOS
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	WB 1:500-1:1000 ICC/IF 1:50-1:200 IP 1:50
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	131117 MW
Uniprot ID	P35228
Expiration Date	12 months from date of receipt.



WB analysis of iNOS expression in Human ...



Immunofluorescent analysis of Raw264.7 c...