



Product Datasheet

Survivin Birc5 Rabbit Monoclonal Antibody (orb547285)


Description	Survivin Birc5 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	A synthesized peptide derived from human Survivin
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	WB 1:500-1:2000IHC 1:50-1:200IP 1:30
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	16298 MW
Uniprot ID	070201
Expiration Date	12 months from date of receipt.



WB analysis of Survivin expression in Neuroblastoma



Immunohistochemical analysis of paraffin-embedded tissue showing Survivin expression in neuroblastoma cells.



Immunohistochemical analysis of paraffin-embedded tissue showing Survivin expression in neuroblastoma cells.