



Product Datasheet

53BP1 TP53BP1 Rabbit Monoclonal Antibody (orb547295)


Description	53BP1 TP53BP1 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, WB
Immunogen	A synthesized peptide derived from human 53BP1
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	WB 1:500-1:2000IHC 1:50-1:200ICC/IF 1:50-1:200FC 1:50
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	213574 MW
Uniprot ID	Q12888
Expiration Date	12 months from date of receipt.



WB analysis of 53BP1 expression in HeLa ...



Immunofluorescent analysis of HepG2 cell...



Immunohistochemical analysis of paraffin...