


Product Datasheet

Integrin alpha 4 ITGA4 Rabbit Monoclonal Antibody (orb547331)

Description	Integrin alpha 4 ITGA4 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthesized peptide derived from human Integrin alpha 4
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	WB 1:500-1:2000
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	114900 MW
Uniprot ID	P13612
Expiration Date	12 months from date of receipt.



WB
analysis of
Integrin
alpha 4
expressi...