

Product Datasheet

Galectin 3 LGALS3 Rabbit Monoclonal Antibody (orb547409)

Description Galectin 3 LGALS3 Rabbit Monoclonal Antibody

Species/Host Rabbit

Reactivity Human, Rat

Conjugation Unconjugated

Tested Applications FC, ICC, IF, WB

Immunogen A synthesized peptide derived from human Galectin 3

Form/Appearance Liquid

Concentration Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Application notes WB 1:500-1:2000 ICC/IF 1:50-1:200 FC 1:50


Isotype Rabbit IgG

Clonality Monoclonal

MW 26152 MW

Uniprot ID [P17931](#)

Expiration Date 12 months from date of receipt.



WB analysis of Galectin 3 expression in ...