

Product Datasheet

p95/NBS1 NBN Rabbit Monoclonal Antibody (orb547458)

Description p95/NBS1 NBN Rabbit Monoclonal Antibody

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications FC, ICC, IF, IHC, IP, WB

Immunogen A synthesized peptide derived from human p95/NBS1

Form/Appearance Liquid

Concentration Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Application notes WB 1:1000-1:2000IHC 1:50-1:200ICC/IF 1:50-1:200IP 1:50FC 1:50


Isotype Rabbit IgG

Clonality Monoclonal


MW 84959 MW

Uniprot ID [O60934](#)

Expiration Date 12 months from date of receipt.



Immunohistochemical analysis of paraffin...



WB analysis of p95/NBS1 expression in He...