

Product Datasheet

MRP2 ABCC2 Rabbit Monoclonal Antibody (orb547554)

Description	MRP2 ABCC2 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, WB
Immunogen	A synthesized peptide derived from human MRP2
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	WB 1:500-1:2000 ICC/IF 1:50-1:200 FC 1:100
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	174207 MW
Uniprot ID	Q92887
Expiration Date	12 months from date of receipt.



WB analysis of MRP2 expression in HepG2 ...