


Product Datasheet

pro Caspase 7 CASP7 Rabbit Monoclonal Antibody (orb547586)

Description	pro Caspase 7 CASP7 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthesized peptide derived from human pro Caspase 7
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	WB 1:500-1:1000IHC 1:50-1:200ICC/IF 1:50-1:200
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	34277 MW
Uniprot ID	P55210
Expiration Date	12 months from date of receipt.



WB analysis of pro Caspase 7 expression

...