



Product Datasheet

AIF/AIFM1 Antibody (monoclonal, 2I5) (orb547790)


Description	AIF/AIFM1 Antibody (monoclonal, 2I5)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, WB
Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human AIF (582-613aa FNRMPiARKIIKDGEQHEDLNEVAKLFNIHED), identical to the related mouse and rat sequences.
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	70 kDa
Uniprot ID	O95831
Expiration Date	12 months from date of receipt.



Western blot analysis of AIF using anti-...



Western blot analysis of AIF using anti-...



IHC analysis of AIF using anti-AIF antib...