

Product Datasheet

CCL19/Mip 3 Beta Rabbit Monoclonal Antibody (orb547805)

Description CCL19/Mip 3 Beta Rabbit Monoclonal Antibody

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Immunogen A synthesized peptide derived from human CCL19

Form/Appearance Liquid

Concentration Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Application notes WB 1:500-1:2000

Isotype Rabbit IgG


Clonality Monoclonal

MW 10993 MW

Uniprot ID [Q99731](#)

Sensitivity 10pg/ml

Expiration Date 12 months from date of receipt.



Western blot analysis of extracts from C...