

## Product Datasheet

# Lysozyme LYZ Rabbit Monoclonal Antibody (orb547865)

**Description** Lysozyme LYZ Rabbit Monoclonal Antibody

**Species/Host** Rabbit

**Reactivity** Human, Mouse, Rat

**Conjugation** Unconjugated

**Tested Applications** ICC, IF, IHC, IP, WB

**Immunogen** A synthesized peptide derived from human Lysozyme

**Form/Appearance** Liquid

**Concentration** Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.

**Storage** Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Note** For research use only

**Application notes** WB 1:5000-1:10000 IHC 1:50-1:200 ICC/IF 1:50-1:200 IP 1:20


**Isotype** Rabbit IgG

**Clonality** Monoclonal


**MW** 16537 MW

**Uniprot ID** [P61626](#)

**Expiration Date** 12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Western blot analysis of Lysozyme using ...