



Product Datasheet

Tyrosine Hydroxylase TH Rabbit Monoclonal Antibody (orb547899)


Description	Tyrosine Hydroxylase TH Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, WB
Immunogen	A synthesized peptide derived from human Tyrosine Hydroxylase
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	WB 1:500-1:2000 ICC/IF 1:50-1:200 FC 1:50
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	58600 MW
Uniprot ID	P07101
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Western blot analysis of Tyrosine Hydrox...



Immunofluorescent analysis of PC-12 cell...