



Product Datasheet

EAAT1 SLC1A3 Rabbit Monoclonal Antibody (orb547958)


Description	EAAT1 SLC1A3 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthesized peptide derived from human EAAT1
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	WB 1:500-1:2000IHC 1:50-1:200
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	59572 MW
Uniprot ID	P43003
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of human br...



Western blot analysis of Mouse brain lys...