

Product Datasheet

ATP citrate lyase ACLY Rabbit Monoclonal Antibody (orb548018)

Description

ATP citrate lyase ACLY Rabbit Monoclonal



Immunohistochemical analysis of paraffin...

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

FC, ICC, IF, IHC, IP, WB

Immunogen

A synthesized peptide derived from human ATP citrate lyase

Form/Appearance

Liquid

Concentration

Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Application notes

WB 1:1000-1:10000 IHC 1:50-1:200 ICC/IF 1:50-1:200 IP 1:50 FC 1:50. Restore with deionized water (or equivalent) for reconstitution volume of 1.0 mL

Isotype

Rabbit IgG

Clonality

Monoclonal

MW


120839 MW

Uniprot ID
[P53396](https://www.uniprot.org/entry/P53396)
Sensitivity

> 5000 cells

Expiration Date

12 months from date of receipt.



Western blot analysis of ATP citrate lya...