



Product Datasheet

NMI Antibody (monoclonal, 2F3) (orb548094)


Description	NMI Antibody (monoclonal, 2F3)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, WB
Immunogen	E.coli-derived human NMI recombinant protein (Position: E2-E307). Human NMI shares 64% amino acid (aa) sequence identity with mouse NMI.
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500µg/ml.
Isotype	Mouse IgG2b
Clonality	Monoclonal
MW	40 kDa
Uniprot ID	Q13287
Expiration Date	12 months from date of receipt.



Western blot analysis of NMI using anti-...



IHC analysis of NMI using anti-NMI antib...



IHC analysis of NMI using anti-NMI antib...