



Product Datasheet

MED15 Antibody (monoclonal, 6F4) (orb548224)


Description	MED15 Antibody (monoclonal, 6F4)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, WB
Immunogen	E. coli-derived human MED15 recombinant protein (Position: M1-A285).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	86 kDa
Uniprot ID	Q96RN5
Expiration Date	12 months from date of receipt.



Western blot analysis of MED15 using anti-MED15 antibody.



IF analysis of MED15 using anti-MED15 antibody.



Flow Cytometry analysis of U2OS cells using anti-MED15 antibody.