


Product Datasheet

EIF2C3 AGO3 Rabbit Monoclonal Antibody (orb548288)

Description	EIF2C3 AGO3 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthesized peptide derived from human EIF2C3
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	WB 1:1000-1:2000
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	97360 MW
Uniprot ID	Q9H9G7
Sensitivity	> 5000 cells
Expiration Date	12 months from date of receipt.



Western blot analysis of EIF2C3 expression...