

Product Datasheet

Caspase 5 CASP5 Rabbit Monoclonal Antibody (orb548357)

Description Caspase 5 CASP5 Rabbit Monoclonal Antibody

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications ICC, IF, IHC, IP, WB

Immunogen A synthesized peptide derived from human Caspase 5

Form/Appearance Liquid

Concentration Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Application notes WB 1:1000-1:2000 IHC 1:50-1:200 ICC/IF 1:50-1:200 IP 1:50


Isotype Rabbit IgG

Clonality Monoclonal


MW 49736 MW

Uniprot ID **P51878**

Expiration Date 12 months from date of receipt.



Western blot analysis of Caspase 5 expression.



IHC analysis of Caspase 5 using anti-Caspase 5 antibody.