

Product Datasheet

beta I Tubulin TUBB1 Rabbit Monoclonal Antibody (orb548367)

Description beta I Tubulin TUBB1 Rabbit Monoclonal Antibody

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications FC, ICC, IF, IHC, WB

Immunogen A synthesized peptide derived from human beta I Tubulin

Form/Appearance Liquid

Concentration Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Application notes WB 1:3000-1:10000 IHC 1:100-1:200 ICC/IF 1:100-1:200 FC 1:100


Isotype Rabbit IgG

Clonality Monoclonal


MW 50327 MW

Uniprot ID [Q9H4B7](#)


Expiration Date 12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Western blot analysis of beta I Tubulin ...



IF analysis of immunocytochemical sectio...