

---

## Product Datasheet

**GM130 GOLGA2 Antibody (monoclonal, 6D4)  
(orb548393)**

**Description**

GM130 GOLGA2 Antibody (monoclonal, 6D4)

**Species/Host**

Mouse

**Reactivity**

Human, Mouse, Rat

**Conjugation**

Unconjugated

**Tested Applications**

FC, ICC, IF, IHC, WB

**Immunogen**

E. coli-derived human GM130 recombinant protein (Position: E796-E913).

**Form/Appearance**

Lyophilized

**Concentration**

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

**Storage**

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

**Note**

For research use only

**Application notes**

Add 0.2ml of distilled water will yield a concentration of 500µg/ml.

**Isotype**

Mouse IgG1

**Clonality**


Monoclonal

**MW**


130 kDa

**Uniprot ID**
[Q08379](https://www.uniprot.org/entry/Q08379)
**Expiration Date**


12 months from date of receipt.



Western blot analysis of GM130 using ant...



Western blot analysis of GM130 using ant...



IHC analysis of GM130 using anti-GM130 a...