

Product Datasheet

Flotillin 2 FLOT2 Rabbit Monoclonal Antibody (orb548402)

Description Flotillin 2 FLOT2 Rabbit Monoclonal Antibody

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications IP, WB

Immunogen A synthesized peptide derived from human Flotillin 2

Form/Appearance Liquid

Concentration Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Application notes WB 1:500-1:2000IP 1:30


Isotype Rabbit IgG

Clonality Monoclonal

MW 47064 MW

Uniprot ID [Q14254](#)

Expiration Date 12 months from date of receipt.



Western blot analysis of Flotillin-2 exp...