



Product Datasheet

Tubulin gamma Rabbit Monoclonal Antibody (orb548412)

Description	Tubulin gamma Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, IP, WB
Immunogen	A synthesized peptide derived from human Tubulin gamma
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	WB 1:500-1:2000IHC 1:50-1:200ICC/IF 1:50-1:200IP 1:50FC 1:50
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	51170 MW
Uniprot ID	P23258
Sensitivity	10pg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Tubulin gamma e...



Immunohistochemical analysis of paraffin...