

Product Datasheet

Histone H2B (mono methyl R79) HIST2H2BE Rabbit Monoclonal Antibody (orb548431)

Description	Histone H2B (mono methyl R79) HIST2H2BE Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthesized peptide derived from human Histone H2B (mono methyl R79)
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	WB 1:1000-1:2000
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	13920 MW
Uniprot ID	Q16778
Expiration Date	12 months from date of receipt.



Western blot analysis of Histone H2B (mo...