

Product Datasheet

Histone H1.0 H1F0 Rabbit Monoclonal Antibody (orb548466)

Description Histone H1.0 H1F0 Rabbit Monoclonal Antibody

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ICC, IF, IHC, WB

Immunogen A synthesized peptide derived from human Histone H1.0

Form/Appearance Liquid

Concentration Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Application notes WB 1:1000-1:2000IHC 1:50-1:200ICC/IF 1:50-1:200

Isotype Rabbit IgG


Clonality Monoclonal

MW 20863 MW

Uniprot ID **P07305**

Sensitivity > 5000 cells

Expiration Date 12 months from date of receipt.



Western blot analysis of Histone H10 exp...