


Product Datasheet

Golgin 97 GOLGA1 Monoclonal Antibody (orb548495)

Description	Golgin 97 GOLGA1 Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthesized peptide derived from human Golgin 97 The Golgi-associated protein golgin A1 (GOLGA1, golgin-97) was first isolated as a Golgi complex autoantigen associated with the autoimmune disorder Sjogren's syndrome. The golgin-97 protein contains a carboxy-terminal GRIP domain and is a commonly used trans-Golgi network (TGN) marker.
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	WB 1:1000-1:5000IP 1:30
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	Q92805
Expiration Date	12 months from date of receipt.



Western blot analysis of GOLGA1 expressi...