


Product Datasheet

His Tag Antibody (monoclonal, 19b10) (orb548525)

| | |
|--------------------------|---|
| Description | His Tag Antibody (monoclonal, 19b10) |
| Species/Host | Mouse |
| Reactivity | All |
| Conjugation | Unconjugated |
| Immunogen | A synthetic peptide corresponding to a sequence of His Tag (HHHHHH). |
| Form/Appearance | Lyophilized |
| Storage | Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles. |
| Note | For research use only |
| Application notes | Add 0.2ml of distilled water will yield a concentration of 500µg/ml. |
| Isotype | Mouse IgG2b |
| Clonality | Monoclonal |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of His Tag using a...