



Product Datasheet

Phospho-Histone H3 (S10) HIST1H3A Rabbit Monoclonal Antibody (orb548537)

Description	Phospho-Histone H3 (S10) HIST1H3A Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, IP, WB
Immunogen	A synthesized peptide derived from human Phospho-Histone H3 (S10).
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	WB 1:1000-1:2000IHC 1:50-1:200ICC/IF 1:50-1:200IP 1:50
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	P68431
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Western blot analysis of Phospho-Histone...