



Product Datasheet

Firefly Luciferase Rabbit Monoclonal Antibody (orb548542)

Description	Firefly Luciferase Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Firefly
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	A synthesized peptide derived from Firefly Luciferase
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	WB 1:500-1:2000 ICC/IF 1:50-1:200
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	60745 MW
Uniprot ID	P08659
Expiration Date	12 months from date of receipt.



Western blot analysis of Firefly Luciferase.



Immunofluorescent analysis of 293 cells ...