

Product Datasheet

Phospho-c-Myc (T58 + S62) Rabbit Monoclonal Antibody (orb548570)

Description Phospho-c-Myc (T58 + S62) Rabbit Monoclonal Antibody

Species/Host Rabbit

Reactivity Human, Rat

Conjugation Unconjugated

Tested Applications FC, ICC, IF, IP, WB

Immunogen A synthesized peptide derived from human Phospho-c-Myc (T58 + S62)

Form/Appearance Liquid

Concentration Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Application notes WB 1:500-1:2000 ICC/IF 1:100-1:200 IP 1:50 FC 1:50


Isotype Rabbit IgG

Clonality Monoclonal

MW 48804 MW

Uniprot ID **P01106**

Expiration Date 12 months from date of receipt.



Western blot analysis of Phospho-c-Myc (...)