

Product Datasheet

Phospho-c-Myc (S62) Rabbit Monoclonal Antibody (orb548571)

Description Phospho-c-Myc (S62) Rabbit Monoclonal Antibody

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ICC, IF, IHC, IP, WB

Immunogen A synthesized peptide derived from human Phospho-c-Myc (S62)

Form/Appearance Liquid

Concentration Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Application notes WB 1:500-1:2000 IHC 1:50-1:200 ICC/IF 100-1:500 IP 1:50


Isotype Rabbit IgG

Clonality Monoclonal


MW 48804 MW

Uniprot ID **P01106**

Expiration Date 12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Western blot analysis of Phospho-c-Myc (...)