


Product Datasheet

Phospho-PDGF Receptor beta (Y740) PDGFRB Rabbit Monoclonal Antibody (orb548580)

Description	Phospho-PDGF Receptor beta (Y740) PDGFRB Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	A synthesized peptide derived from human Phospho-PDGF Receptor beta (Y740)
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	WB 1:1000-1:2000 ICC/IF 1:50-1:200
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	123968 MW
Uniprot ID	P09619
Expiration Date	12 months from date of receipt.



Western blot analysis of Phospho-PDGF Re...