



Product Datasheet

Phospho-Tau (S396) MAPT Rabbit Monoclonal Antibody (orb548582)

Description	Phospho-Tau (S396) MAPT Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthesized peptide derived from human Phospho-Tau (S396)
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	WB 1:10000-1:20000IHC 1:100-1:500
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	78928 MW
Uniprot ID	P10636
Expiration Date	12 months from date of receipt.



Western blot analysis of Phospho-Tau (S396)



Immunohistochemical analysis of paraffin...