


Product Datasheet

Phospho-YAP1 (S127) Rabbit Monoclonal Antibody (orb548585)

Description	Phospho-YAP1 (S127) Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthesized peptide derived from human YAP1
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	WB 1:5000-1:20000IHC 1:50-1:200
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	54462 MW
Uniprot ID	P46937
Expiration Date	12 months from date of receipt.



Western blot analysis of Phospho-YAP1 (S...

