



Product Datasheet

Phospho-TrkB (Y817) NTRK2 Rabbit Monoclonal Antibody (orb548631)

Description	Phospho-TrkB (Y817) NTRK2 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, IP, WB
Immunogen	A synthesized peptide derived from human Phospho-TrkB (Y817)
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	WB 1:500-1:1000IHC 1:50-1:100ICC/IF 1:50-1:100IP 1:30
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	91999 MW
Uniprot ID	Q16620
Expiration Date	12 months from date of receipt.



Western blot analysis of Phospho-TrkB (Y817)



Immunohistochemical analysis of paraffin-embedded tissue sections