



Product Datasheet

Phospho-Smad5 (S463/465) Rabbit Monoclonal Antibody (orb548636)

Description	Phospho-Smad5 (S463/465) Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthesized peptide derived from human Phospho-Smad5 (S463/465)
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	WB 1:500-1:2000IHC 1:50-1:200ICC/IF 1:50-1:200
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	52258 MW
Uniprot ID	Q99717
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Western blot analysis of Phospho-Smad5 i...