

Product Datasheet

Phospho-Histone H1.4 (T17) HIST1H1E Rabbit Monoclonal Antibody (orb548656)

Description Phospho-Histone H1.4 (T17) HIST1H1E Rabbit Monoclonal Antibody

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ICC, IF, IHC, WB

Immunogen A synthesized peptide derived from human Phospho-Histone H1.4 (T17)

Form/Appearance Liquid

Concentration Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Application notes WB 1:500-1:2000 IHC 1:50-1:200 ICC/IF 1:50-1:200


Isotype Rabbit IgG

Clonality Monoclonal


MW 21865 MW

Uniprot ID **P10412**, **P16402**

Expiration Date 12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Western blot analysis of Phospho-Histone...