


Product Datasheet

Phospho-Histone H2B (T129) HIST1H2BC Rabbit Monoclonal Antibody (orb548663)

Description	Phospho-Histone H2B (T129) HIST1H2BC Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Yeast
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthesized peptide derived from yeast Phospho-Histone H2B (T129)
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	WB 1:500-1:1000
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	14252 MW
Uniprot ID	P02293
Expiration Date	12 months from date of receipt.



Western blot analysis of Histone H2B exp...