

## Product Datasheet

### Bovine IFN $\alpha$ ELISA Kit (Ready to Use) (orb548690)

<b>Description</b>	Bovine Interferon Alpha (IFN $\alpha$ ) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for co...
<b>Reactivity</b>	Bovine
<b>Range</b>	15.6-1000pg/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q04790</a>
<b>Sensitivity</b>	5.9pg/mL
<b>Entrez</b>	<a href="#">282257</a>
<b>Expiration Date</b>	Please enquire.