

Product Datasheet

Bovine ITaC ELISA Kit (Ready to Use) (orb548693)

| | |
|--------------------------|---|
| Description | Bovine Interferon Inducible T-Cell Alpha Chemoattractant (ITaC) Sandwich ELISA KitReady to Use: Pre-... |
| Reactivity | Bovine |
| Range | 46.8-3000pg/mL |
| Form/Appearance | Ready to Use |
| Storage | All reagents should be stored at -20°C upon receipt |
| Note | For research use only |
| Application notes | Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids. |
| Sample Types | serum, plasma and other biological fluids. |
| Uniprot ID | A9QWQ1 |
| Sensitivity | 15.46pg/mL |
| Entrez | 516104 |
| Expiration Date | Please enquire. |