

Product Datasheet

Bovine IL13 ELISA Kit (Ready to Use) (orb548694)

Description	Bovine Interleukin 13 (IL13) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conv...
Reactivity	Bovine
Range	31.2-2000pg/ml
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q9XSV9
Sensitivity	13.4pg/ml
Entrez	281247
Expiration Date	Please enquire.