

Product Datasheet

Canine ALB ELISA Kit (Ready to Use) (orb548704)

Description	Canine Albumin (ALB) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient a...
Reactivity	Canine
Range	1.56-100µg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P49822
Sensitivity	0.49µg/mL
Entrez	403550
Expiration Date	Please enquire.