

Product Datasheet

Canine IGF1 ELISA Kit (Ready to Use) (orb548705)

Description	Canine Insulin Like Growth Factor 1 (IGF1) Sandwich ELISA KitReady to Use: Pre-diluted detection rea...
Reactivity	Canine
Range	2.5-160ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P33712
Sensitivity	0.9ng/mL
Entrez	610255
Expiration Date	Please enquire.